

## **LIFEPORT CMS the choice for Air New Zealand**

The LifePort Commercial Medical Stretcher (CMS) has been delivered and approved by FAA STC for the full range of Air New Zealand's fleet, including the Airbus A320, Boeing 747, 767, and 777-200ER.

“For our operation, there really was no better choice than the LifePort stretcher,” said Air New Zealand Engineering Services General Manager Chris Nassenstein. “We operate a fleet of mixed aircraft types, and the CMS meets all these requirements, as well as being extremely easy to install and adjust. In addition, the CMS bed top is wider and longer than anything else on the market, providing a higher level of comfort than other stretchers”.

Some highlights of the LifePort Commercial Medical Stretcher include:



- 78” long, 30” wide bed top with a plush 3” thick mattress
- Curved support legs to contour around seats in restricted pitch settings
- Adjustable leg attachment points, allowing installation across full spectrum of commercial aircraft
- Height-adjusted privacy curtain assembly to fit range of aircraft overhead bins
- Medical transport system, with a variety of customer-specific accessories, including incubator transport option
- Built to TSO, certified by FAA STC
- Assembled in sections to facilitate installation, removal, and storage

“LifePort is thrilled to deliver these stretchers to Air New Zealand,” says Gary Darley, CEO of LifePort, Inc. “Air New Zealand has been a valuable resource for us in the course of this project, providing product feedback and industry guidance which will exponentially improve LifePort’s ability to service similar customers. Having an airline like Air New Zealand select the CMS after a direct comparison with our main competitor re-enforces our belief that the CMS is the best stretcher on the market. We are very proud of this product and its ability to fit into existing airline operations – it is always gratifying when a customer and operator share that opinion.”

For more information contact:

**Sales Department**  
**1610 Heritage Street**  
**Woodland, WA 98674 USA**  
**Tel: (360) 225-1212 Fax: (360) 225-1214**  
[www.lifeport.com](http://www.lifeport.com)  
[solutions@lifeport.com](mailto:solutions@lifeport.com)